

REXMED™

試管振盪器

VORTEX MIXER

MEZCLADOR DE VÓRTICE

MÉLANGEUR DE VORTEX

VORTEX-MISCHER

ВИХРЕВОЙ СМЕСИТЕЛЬ



ISO 9001
ISO 13485



RVM-102 **Designed For Better Life**

APPLICATIONS

--This equipment is designed for mixing solutions in test tubes or small flasks by means of a circular, oscillating movement.

FEATURES

- Compact, robust and durable designed for efficiency.
- Variable speed control.
- Adopting high quality motor to achieve the performance.
- Continuous operation or automatically press star.
- 3 rubber feet for quite vortex mixing action and high stability of the unit avoiding vibrations transmission to the lab bench.
- Rubber platform head is designed for easy cleaning and maintenance.

CONTROL AND OPERATION

- Placing this unit on a firm level surface and connect to the main power supply.
- Load samples to be analyzed and mixed on the top head.
- Simply and easy switch to AUTO for continuous operation or switch to TOUCH for automatically press star function.
- Samples can be load or unload without turning off the unit.
- Speed can be adjusted from 0~3,000rpm.



SPECIFICATIONS

| | | | |
|---------------|-------------|---------------|-----------------|
| Model | RVM-102 | Top head type | Rubber platform |
| Voltage | AC 110/220V | N.W. | 4.5 KGS |
| Speed | 3,000 rpm | G.W. | 5 KGS |
| Top head size | Ø10 cm | Packing size | 18 x 14 x 18 cm |



* Specifications subject to change without prior notice due to product improvement.
 * Manufactured by ISO 9001, ISO 13485, GMP registered company.

Copyright © since 1976 REXMED INDUSTRIES CO., LTD.
 TEL: 886-7-3814010 FAX: 886-7-3814013
 E-MAIL: sales@rexmed.com WEB: www.rexmed.com, www.rexmed.eu