



COLONY COUNTER - GMP Model

Applications:

- ❖ **LAPIZ** Digital Colony Counter is ideal for fast and fatigue free counting of Bacterial Colonies.

LAPIZ Advantages:

- ❖ Circular tube and powerful magnifying lens - reduces operator fatigue.
- ❖ Streamlined design ensures fast and accurate count.
- ❖ Desk top compact model.

Construction:

- ❖ Re-designed model with advanced sensors enables usage of any commercially available hard tipped marker pen.
- ❖ LED Display with audio beep at every count.
- ❖ Wouffhuegel grid glass plate.
- ❖ On marking the colony by applying a slight pressure with the marker pen, a count is registered and confirmed by an audio beep.
- ❖ Full Stainless steel fabricated body duly heat cured epoxy coated from outside for better aesthetics.



Features:

- ❖ Magnification: 2x.
- ❖ Maximum Count: 9,999.
- ❖ Manual count switch for counting without "Marker Pen".

Model	Product	Maximum Count	App. Shipping Details Net / Gr. / Volume
#0671M	Digital - Desk Model - SS GMP Model	9,999	4.5 kg / 5.0 kg (410 x 320 x 300 mm) 0.04 m ³

All electrical appliances work on 230V AC, 50Hz, single phase. Other Voltages available on request.
Due to continuous improvement in design, the product supplied may have modified features.

Our Other Products: Autoclave Vertical: Standard / PAD / Fully Automatic / Jacketed, Autoclave Portable: Standard / Automatic, Autoclave SLL Series: Vacuum, Autoclave SLE series • Media Preparator • Table Top Autoclaves: Normal / Class B • UV Sterilization Cabinet • Rotary Vacuum Evaporator • Vacuum Controller • Chilled Water Circulator • PTFE Diaphragm Vacuum Pump • BOD Incubator • BSC Clean Bench • Water Still / Stone-fin / Automatic Water Still • Water Bath: Stirred (Constant) / Serological / Concentric Rings / Refrigerated Circulating Bath • Round Bath: Clear Bottom / Immersion Series / ECO Series • Ultrasonic Bath: Digital / LCD • Dry Bath • Steam pot • Colony Counter • Hot Plate with Magnetic Stirrer • Anaerobic Jar: SS / Polycarbonate • Pipette & Petri Dish Sterilizing Box • Test Tube Rack • ESR Westergren & Wintrobe Apparatus • Slide Staining Stand / Rack / Trough • Slide Trays • Slide Cabinet



MEDICA INSTRUMENT MFG. CO.

105, Vasan Udyog Bhavan, Sun Mill Compound, Lower Parel (West)
Mumbai 400013, INDIA • ☎ +91-22/6618-9999 • 📠 +91-22/6618-9996
✉ sales@medicainstrument.com • www.medicainstrument.com

