



For reliable steaming / evaporation

EQUITRON®

CONCENTRIC RINGS WATER BATH

Digital

Application:

- ❖ EQUITRON Concentric Rings Water Baths are ideal for all applications requiring gentle steaming work with Erlenmeyer flasks, glasses etc.

EQUITRON Advantages:

- ❖ Fitted with a water level indicating device.
- ❖ Optional constant water level device.
- ❖ A removable flat drop on lid for concentric rings enables the bath to be used for other work also.

Construction:

- ❖ Double walled, insulated, full stainless steel construction.
- ❖ Internal chamber, external chamber, drop on lid with holes of 75mm Ø all of stainless steel.
- ❖ Concentric rings in steps of 25 and 50 mm ID moulded from industrial grade non-corrosive material.

Standard Features:

- ❖ A precise Digital Temperature Indicator cum Controller having resolution and control accuracy of $\pm 1^\circ\text{C}$ in digital model.
- ❖ Temperature range: Ambient $+5^\circ\text{C}$ to 100°C .
- ❖ Industrial grade long life special heaters for efficient transfer of heat.



Model	Bath	Control	Concentric Rings	Heater	App. Shipping Details Net / Gr. / Volume
#6806DI	30 x 20 x 10 cm	Digital	6 holes of 75 mm Ø with steps of 25 mm	800 W x 1	5.0 kg / 5.6 kg / 0.04 m ³
#6812DI	40 x 30 x 10 cm	Digital	12 holes of 75 mm Ø with steps of 25mm	800 W x 2	9.0 kg / 10 kg / 0.06 m ³
#6824DI	60 x 40 x 10 cm	Digital	24 holes of 75 mm Ø with steps of 25mm	1500 W x 2	13 kg / 28 kg / 0.15 m ³

All electrical appliances work on 230 VAC, 50Hz, single phase.
Due to continuous improvement in design, the product supplied may have modified features.

Our Other Products: Autoclave Vertical: Standard / PAD / Fully Automatic / Jacketed, Autoclave Portable: Standard / Automatic, Autoclave SLL Series: Vacuum, Autoclave SLE series • Media Preparator • Table Top Autoclaves: Normal / Class B • UV Sterilization Cabinet • Rotary Vacuum Evaporator • Vacuum Controller • Chilled Water Circulator • PTFE Diaphragm Vacuum Pump • BOD Incubator • BSC Clean Bench • Water Still / Stone-fin / Automatic Water Still • Water Bath: Stirred (Constant) / Serological / Concentric Rings / Refrigerated Circulating Bath • Round Bath: Clear Bottom / Immersion Series / ECO Series • Ultrasonic Bath: Digital / LCD • Dry Bath • Steam pot • Colony Counter • Hot Plate with Magnetic Stirrer • Anaerobic Jar: SS / Polycarbonate • Pipette & Petri Dish Sterilizing Box • Test Tube Rack • ESR Westergren & Wintrobe Apparatus • Slide Staining Stand / Rack / Trough • Slide Trays • Slide Cabinet



MEDICA INSTRUMENT MFG. CO.

105, Vasan Udyog Bhavan, Sun Mill Compound, Lower Parel (West)
Mumbai 400013, INDIA • ☎ +91-22/6618-9999 • 📠 +91-22/6618-9996
✉ sales@medicainstrument.com • www.medicainstrument.com



0010_1005-3-1