

Various positions of RAU-OT 8005

REMOTE CONTROL ELECTRO O.T. TABLE WITH LONGITUDINAL TOP SLIDE



RAU-OT 8005

TECHNICAL DATA

Top Slide both Head & leg side: 30 cm
 Minimum Height: 75 cm
 Maximum Height: 105 cm Length of Table Top: 215 cm
 Width of Table Top: 50 cm
 Lateral Tilt: + 20° / -20°
 Trendelenburg: +25° / -25°
 Back Section: +90° / -40°
 Head & leg Section: +90° / -90° manual
 Operated Kidney Elevator: 15 cm manual

SALIENT FEATURES:

100% Life Saving due to low voltage Electricity Supply.
 Offline Battery Back-Up
 360kg (800 lb) weight capacity
 Full Table Functions
 Full Articulation Capability
 Removable Leg section
 Full Body imaging without special accessories
 Durable one-touch pendant control
 Electrically Controlled Height
 Trendelenburg/Reverse, Lateral Tilt.
 4-Section seamless moulded mattress
 Support on a heavy base column covered with stainless steel sheet.

ACCESSORIES OF O.T. TABLE

Anaesthetic Frame (1Pcs)
 Padded Side Supports (1 Pair)
 Intravenous Armboard (1 Pair) Padded
 Shoulder Supports (1 Pair) Knee
 Crutches (1 Pair)
 One Mattress Set Thickness 50mm

