



SURGICAL & EXAMINATION LAMP



ISO 9001

ISO 13485

ROL-101S **Designed For Better Life**



REXMED Industries Co., Ltd.

P.O.BOX 487, Kaohsiung, Taiwan

Tel: +886-7-3814010 Fax: +886-7-3814013

E-mail: info@rexmed.com, sales@rexmed.com

<http://www.rexmed.com>



ROL-101S



- Main construction: powder coated steel and ABS.
- Mobile stand type 5 castors with 3 lockable castors.
- The halogen lamp adopts snake hose type that can make directions and angles be freely adjusted to fulfil multi-functional illumination.
- The halogen lamp can gather light and form in about 12~20cm diameter light circle.
- Stand-type designed, light weight and easy installation. The halogen lamp is suitable for the emergency ward, dispensary, dressing room, examination room, etc.
- Its inside is installed with one transformer that can transform AC 110/230V into DC 12V for safety use.

SPECIFICATIONS

Power	AC 110/230V	Controller	No/Off switch
Lamp head	7cm	Snake hose	77cm
Bulb No.	OSRAM 12V 35W x 1	Lamp head elevation	60 cm
Bulb life time	4,000 hrs \pm 10%	N.W.	5 KGS
Intensity	13,000 Lux/1m	G.W.	6 KGS
Colour Temperature	3,500K \pm 200K	Packing size	L60 X W60 X H16cm