

Ultrasonic, Ultra-Clean

Sturdy Ultrasonic Cleaners can burst waves up to 46,000 cycles per second, producing longitudinal wave forms that cause a build up of dielectric pressure, which results in micro-vacuum bubbles that explode with forceful cleansing action, easily removing dirt and grime from the most hard-to-reach corners and crevices.



| Model | Internal Tank Size (LxWxH) | External Dimensions | Capacity | Tank | Cover | Basket | Housing Material | Ultrasonic Output | Heater Output | Electronic Timer |
|--------|----------------------------|---------------------|----------|--------------|--------------|--------------------------|---------------------------------|-------------------|---------------|------------------|
| UC-400 | 495x295x150 mm | 530x325x400 mm | 22L | St. St. #304 | St. St. #304 | St.St. #304 | Aluminum alloy (powder coating) | 400W / 40 KHz | 400W | 1-99 mins. |
| UC-300 | 295x240x150 mm | 320x265x280 mm | 10.5L | St. St. #304 | St. St. #304 | St.St. #304 | Aluminum alloy (powder coating) | 300W / 46 KHz | 400W | 1-30 mins. |
| UC-150 | 295x150x100 mm | 320x175x240 mm | 4.5L | St. St. #304 | St. St. #304 | St.St. #304 | Aluminum alloy (powder coating) | 150W / 46 KHz | 200W | 1-30 mins. |
| UC-80 | 150x135x100 mm | 175x160x210 mm | 2L | St. St. #304 | St. St. #304 | St.St. #304 | Aluminum alloy (powder coating) | 80W / 46 KHz | N/A | 1-30 mins. |
| UC-50 | 168x89x55 mm | 193x116x118 mm | 0.6L | St. St. #304 | PS Plastic | ABS Plastic ² | ABS Plastic ² | 50W / 46 KHz | N/A | N/A |

² U/L Fire Resistant, Power Switches on all machines are soft touch buttons, Power Source 110V or 220V or 240V (as requested on order)