



No Heater Burnouts

EQUITRON[®]

WATER STILL

With Thermostatic cut-off

Applications:

- ❖ **EQUITRON** Water Still is ideal for all laboratories where distill water, free from dissolved salts, clear, odourless and tasteless water is required.

EQUITRON Advantages:

- ❖ Thermostatic cut-off (self resetting type) for assured low water detection - with visual indication.
- ❖ Unique water level maintaining device improves efficiency.
- ❖ Compact & user friendly space saving design - occupies less space & can be table or wall mounted.
- ❖ Easily openable Boiling Chamber as well as Condenser for easy cleaning.

Construction:

- ❖ Fully stainless steel.
- ❖ Press fit boiling chamber with internal gasket & condenser.



Standard Accessories:

- ❖ Silicone tubing for distill water output + for flow of water from Condenser to Boiling Chamber.
- ❖ PVC tubing for waste water + raw water with tap adapter and a hose clip.
- ❖ Flow rate meter plus Hoffman clip for raw water flow monitoring.
- ❖ Mounting screws for Water Still and Flow meter.

Model	Capacity	Heater	App. Shipping Details Net / Gr. / Volume
#7655M	4 lit / hr	3.0kW	4.5 kg / 5.5 kg / 0.069 m ³

All electrical appliances work on 230V AC, 50Hz, single phase. Other Voltages available on request.
Due to continuous improvement in design, the product supplied may have modified features.

Our Other Products: Autoclave Vertical: Standard / PAD / Fully Automatic, Autoclave Portable: Standard / Automatic, Autoclave SLL Series: Vacuum, Autoclave SLE series
• Media Preparator • Table Top Autoclaves: Normal / Class B • UV Sterilization Cabinet • Rotary Vacuum Evaporator • Vacuum Controller • Chilled Water Circulator • PTFE Diaphragm Vacuum Pump • BOD Incubator • BSC Clean Bench • Water Still / Stone-fin / Automatic Water Still • Water Bath: Stirred (Constant) / Serological / Concentric Rings / Refrigerated Circulating Bath • Round Bath: Clear Bottom / Immersion Series • Ultrasonic Bath: Digital / LCD • Dry Bath • Steam pot • Colony Counter • Hot Plate with Magnetic Stirrer • Anaerobic Jar: SS / Polycarbonate • Pipette & Petri Dish Sterilizing Box • Test Tube Rack • ESR Westergren & Wintrobe Apparatus • Slide Staining Stand / Rack / Trough • Slide Trays • Slide Cabinet



MEDICA INSTRUMENT MFG. CO.

105, Vasan Udyog Bhavan, Sun Mill Compound, Lower Parel (West)
Mumbai 400013, INDIA • ☎ +91-22/6618-9999 • 📠 +91-22/6618-9996
✉ sales@medicainstrument.com • www.medicainstrument.com



0015_0511-2-1